


# DDC-8/BlackBox QuickStart User's Guide

- 1) Turn on the DDC-8 rack module. (You might hear **3 ding sounds** from the computer speakers)
- 2) Click on the **BlackBox ANSEL** icon to open the DAQ software.  
Then choose **DDCT>Main** from the top menu bar.
- 3) Begin the initialization procedure and
  - a. click on Initialize USB (But **only once!**). The **Log** window should display "...**Interface: 0**"
    - i. If you see something else (possibly -5 or -22), close Blackbox, push the button labeled "**Reset USB**" on the front panel of the DDC-8 module (or power-cycle the device)
  - b. Then choose **DDCT>SDAC** from the top menu and set signal amplitude offsets for each channel in the boxes on the right-hand side. Click on button "**Set SDAC**" to apply the settings.
    - i. In most spectroscopy measurements, the number "**32**" should be entered for each of the analog inputs to be used. This enables full scale digitization of ~2V pulses.
    - ii. In principle, other offsets can be used if specifically required (e.g. if negative [use "**225**"] or bipolar [use "**128**"] pulses are expected.).
- 4) Recheck on the scope. the **form** of the positive **analog signal** (→**Ch 1,..**) and its timing rel. the negative digital trigger pulse (→**NIM\_IN\_0** BNC, no 50-Ω terminator)
- 5) Go back to the **DDCT>Main** window and set following values in the dialog boxes:
  - a. Check the box for "**capture peaks**", to measure amplitudes.
    - i. For initial setup test, choose a **value between 2-5** for the **Number of Events per Buffer**. (This determines how often the GUI display will update information. It can be increased in actual data taking runs.)
  - b. Set each of the desired Ch#s to **max** and **ena**.
  - c. Press **Apply** to save all of the settings adjusted so far.
  - d. Adjust the **Post-trigger delay** to **80 for initial setup**.
    - i. This instructs the DDC-8 about the number of signal samples (40ns ea.) to inspect. A good start value is **80**, but typical larger values are used, like **512**.
  - e. Press the button "**Send delay to board**" to save the delay.
- 6) Answer the "**Get N events**" question box, choose "**0**" (zero) for indefinite run time.
- 7) Press the "**Start**" button to initiate the DDC-8. The counters "**Run time**" and "**Events captured**" should start incrementing, indicating proper functioning of the DAQ.
- 8) The DDC-8 DAQ provides an option to **monitor the event-by-event data stream**, without disturbance, by displaying a corresponding frequency (pulse-height) spectrum as a histogram on an extra window with manipulation options. Proceed as follows:
  - a. From the top menu, select **DDCT>Histograms** to display the window **histograms**, which provides access and manipulation to up to 8 histograms (Hist0,...,Hist7).
  - b. Configure the histogram parameters **Number of Bins** to display, number of **Offset** ADC-channels, **Bin** size in ADC channels, which define the appearance of the histogram in the display
    - i. Default settings are **4096** Bins, **2048** channels offset & **16** channels per bin. The screenshot below has 4096, 0, and 8 instead for displaying data entering DDC8-Ch1.
    - ii. Comparisons are possible by displaying the same data on different histograms, e.g., on Hist0 and on Hist3.
  - c. Initialize the histogram by pressing the button **Init Hist0,..**
  - d. Histogram display is useful for monitoring data taking runs. Event data remain intact.
  - e. Useful hint: Identify zero pulse-height by briefly disconnecting the analog signal cable from the corresponding DDC8 input connector, e.g., Ch1,..
- 9) Once the above data monitoring suggests proper DAQ settings, event data can be written to file by checking the "**Write Peaks to File**" box on the **Main** window (below Start and Stop buttons).



## Tips/Troubleshooting:

- 1) After changing any setting, you need to press "**Apply**", then "**Send delay to board**" before starting data collection. This finalizes the changes in the software, then sends them to the FPGA.
- 2) Be sure to write down in your logbook the Runtime & # of events collected after each run.
  - a. Don't worry if you forget, however, as this information is accessible from the output file.
- 3) Be sure to save all of your data on the EXPERIMENTAL\_DATA partitions on the hard drives (Disk *E*). Most of the C: drives will fill up quickly.
  - a. This was a common problem in previous years. If your data is not saved correctly, be sure to check and make sure the HDD isn't full.
- 4) Keep an eye on the Log Window. If the machine suddenly starts to display "-5"s, there is likely a problem with the connection and it is better to power cycle the instrument and reinitialize.
- 5) Sometimes, you will hear the "ding" of the USB disconnect. This happens occasionally via static discharge or jostling of cables. Power cycle the box & reinitialize to continue
- 6) If you want to quickly export the histogram from Blackbox, double click on the legend entry for "**Hist\_0**" and a window will appear with the bin values in a list. This makes for easy importation into IGOR, but you should practice making the histograms manually on your own.
- 7) If you want to adjust the scaling on the histogram to zoom to a particular region, double click on either axis to bring up the menu.
  - a. Be sure to uncheck the "**autoscale**" box in order to change the scale of a particular axis.